INSTALLATION INSTRUCTIONS

These instructions with the specific connecting diagram delivered with this product are designed to ensure its correct integration into the specified dental unit. They must be followed to meet all applicable requirements. The product must only be built in the dental unit by qualified technicians authorized by its manufacturer. It is the installer's responsibility to design and assemble the dental unit according to applicable standards, to conduct final testing in the appropriate environment, and to deliver information how to use the product as a feature of the dental unit.

INFORMATION

The P18 is an ultrasonic generator for dental scaling fully compatible with the EMS Piezon® handpieces.

The P18 is homologated according to the standards for electrical security EN 60601-1 including the deviations for United States and Canada (CB-Certificates EN 60601-1, ANSI/AAMI ES STD 60601-1, CAN/CSA STD C22.2 NO.60601-1).

According to the electromagnetic conformity standards EN 60601-1-2 this ultrasonic generator is guaranteed not to disturb the EMC of the unit. The tests were passed without need of additional component.

Use of an EN 60601-1 conform power supply as well as the use of EMS approved handpieces and cords are required to retain the conformity with this standard.

The P18 can be built-in dental units without further investigation about the generator itself.

LAYOUT

Install the ultrasonic generator into the dental unit using the hardware provided. Place the ultrasonic generator in a ventilated area and away from any heat source. The ultrasonic generator must be installed in a dry zone as its housing does not provide protection against water or liquids.

CONNECTIONS

Please follow the specific connecting diagram delivered by EMS for the specified dental unit.

Handpiece cord

It is recommended to provide a connector on the dental unit for the handpiece cord. Install a direct connection between the connector and the dental unit water inlet. The anti-extension string must be securely tied up to the connector, on unit side, or to the dental unit if a connector is not used.

SAFETY PRECAUTIONS

EMS accepts no liability for direct or consequential injury or damage resulting from improper installation or use arising in particular through non-observance of EMS instructions or improper preparation and maintenance.

 This product must only be installed and used by trained and qualified personnel. It must exclusively be repaired by an approved EMS repair center.

 Do not attempt to repair or modify this product without an authorization from the manufacturer.

 Contraindications: ultrasonic oscillations may prevent cardiac pacemakers and defibrillators functioning properly. EMS recommends that patients with a cardiac pacemaker should not be treated with this product.

 Risk of explosion: Do not use this product in the presence of anesthetics or inflammable gas.

 Risk of electric shock for the patient or the user: This product must only be connected to the correct power supply, otherwise its safety and correct operation cannot be guaranteed.

 This product has been investigated with regard to safety from electrical shock and fire hazard. ETL (Electrical Testing Laboratories) has not investigated other physiological effects.

 The warranty of this product will be cancelled if you try to open it. The screws on the generator housing must under no circumstance be tampered with.

About this manual

Note that the English version of this manual is the master from which translations derive. In case of any discrepancy, the binding version is the English text.

Technical modifications

EMS reserves the right to modify the technique, accessories, operation instructions or contents of the set due to technical or scientific improvements.

 Typical illustrations, the "Packing list" shows exactly what the product includes depending on the application.

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