

SmartPeg Reference Guide

Valid from September 2017 | 25016-65

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New	
A.B. DENTAL DEVICES	Screw Type, Conical Groovy				3.75		100440	Type 32	
	Screw Type, Conical Groovy				4.5		100431	Type 27	
ADDENTA	Bifasico				4.0		100374	Type 14	
ADIN DENTAL IMPLANT SYSTEM	One				3.3		100390	Type 17	
	One				4.2		100390	Type 17	
	Swell				3.3		100480	Type 49	
	Swell				3.75		100480	Type 49	
	Swell				4.2		100480	Type 49	
	Swell				5.0		100480	Type 49	
	Swell				6.0		100480	Type 49	
	Touareg CloseFit				NP		100545	Type 62	
	Touareg CloseFit				RP		100543	Type 60	
	Touareg CloseFit				WP		100480	Type 49	
	Touareg-S				4.2		100480	Type 49	
Touareg-S				5.0		100480	Type 49		
AGS MEDIKAL	Implance Dental Implant System, Bone Level						100378	Type 06	
	Implance Dental Implant System, Mini Type						100455	Type 38	
	Implance Dental Implant System, Tissue Level						100378	Type 06	
ALPHA BIO	ATI						100425	Type 26	
	ATID				3.75		100425	Type 26	
	Neo CS				3.75		100583	Type 78	New
	NICE 3.2 Conical connection				3.2		100554	Type 63	
	SPI						100425	Type 26	
ALPHA-DENT	Active 3D				3.25		100480	Type 49	
	Active 3D				3.75		100480	Type 49	
	Active 3D				4.2		100480	Type 49	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
ALPHA-DENT	Active 3D				5.0		100480	Type 49
	Active Conus				3.25		100480	Type 49
	Active Conus				3.75		100480	Type 49
	Active Conus				4.2		100480	Type 49
	Active Conus				5.0		100480	Type 49
	Active Plus				4.2		100480	Type 49
	Active Plus				5.0		100480	Type 49
	Classic 3D				3.3		100480	Type 49
	Classic 3D				3.75		100480	Type 49
	Classic 3D				4.2		100480	Type 49
	Classic 3D				5.0		100480	Type 49
	Classic Conus				3.3		100480	Type 49
	Classic Conus				3.75		100480	Type 49
	Classic Conus				4.2		100480	Type 49
	Classic Conus				5.0		100480	Type 49
ALPHA-GATE	Bioactive				3.75		100446	Type 35
ALPHATECH					3.4		100477	Type 46
					3.8		100477	Type 46
					4.3		100477	Type 46
					5.0		100477	Type 46
					6.0		100477	Type 46
					7.0		100477	Type 46
ANKYLOS	C/X		A				100388	Type 16
	C/X		B				100388	Type 16
	C/X		C				100388	Type 16
ANTHOGYR	Anthofit External Hexagon HE						100353	Type 01
	Anthofit Internal Octagon OI						100398	Type 18
	Axiom PX				3.4		100455	Type 38
	Axiom PX				4.0		100455	Type 38
	Axiom PX				4.6		100455	Type 38
	Axiom PX				5.2		100455	Type 38

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
ANTHOGYR	Axiom REG			3.4		100455	Type 38	
	Axiom REG			4.0		100455	Type 38	
	Axiom REG			4.6		100455	Type 38	
	Axiom REG			5.2		100455	Type 38	
	Ossfit					100350	Type 04	
ARDS IMPLANTS	Classic			All		100431	Type 27	
	Premium			All		100431	Type 27	
	Smart			All		100431	Type 27	
AVINENT IMPLANT SYSTEM	Coral External			3.5		100436	Type 30	
	Coral External			4.1		100353	Type 01	
	Coral External			5.1		100357	Type 03	
	Internal Hexagon Coral			3.5		100440	Type 32	
	Internal Hexagon Coral			4.1		100440	Type 32	
	Internal Hexagon Ocean			3.5		100440	Type 32	
	Internal Hexagon Ocean			4.1		100440	Type 32	
	Ocean External			3.5		100436	Type 30	
	Ocean External			4.1		100353	Type 01	
AXIS	Ceramic			4.0		100455	Type 38	
B&B DENTAL	Duravit SLIM Implant			3.0		100404	Type 22	New
	Duravit SLIM Implant			3.4		100404	Type 22	New
	Duravit 3P Implant			3.5		100425	Type 26	New
	Duravit 3P Implant			4.0		100425	Type 26	New
	Duravit 3P Implant			4.5		100425	Type 26	New
	Duravit 3P Implant			5.0		100425	Type 26	New
	DURAVIT EV Implant			4.5		100425	Type 26	New
	DURAVIT EV Implant			4.0		100425	Type 26	New
	DURAVIT EV Implant			4.5		100425	Type 26	New
	DURAVIT EV Implant			5.0		100425	Type 26	New
	DURAVIT WIDE Implant			5.0		100425	Type 26	New
	DURAVIT WIDE Implant			6.0		100425	Type 26	New
	DURAVIT WIDE Implant			6.5		100425	Type 26	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
B&B DENTAL	DURAVIT WIDE Implant			7.0		100425	Type 26	New
BEGO IMPLANT SYSTEMS	Semados Mini			2.7		100355	Type 02	
	Semados Mini			2.9		100355	Type 02	
	Semados Mini			3.1		100355	Type 02	
	Semados RI-line			3.75		100425	Type 26	
	Semados RI-line			4.1		100425	Type 26	
	Semados RI-line			4.5		100425	Type 26	
	Semados RI-line			5.5		100425	Type 26	
	Semados RS/RSX			3.0		100481	Type 50	
	Semados RS/RSX			3.75		100425	Type 26	
	Semados RS/RSX			4.1		100425	Type 26	
	Semados RS/RSX			4.5		100425	Type 26	
	Semados RS/RSX			5.5		100425	Type 26	
	Semados S-line			3.25		100425	Type 26	
	Semados S-line			3.75		100425	Type 26	
	Semados S-line			4.1		100425	Type 26	
	Semados S-line			4.5		100425	Type 26	
	Semados S-line			5.5		100425	Type 26	
BIO3	Implant Advanced ATI			3.3		100480	Type 49	
	Implant Advanced ATI			3.8		100480	Type 49	
	Implant Advanced ATI			4.2		100480	Type 49	
	Implant Advanced ATI			5.0		100480	Type 49	
	Implant Progressive PTI			3.3		100480	Type 49	
	Implant Progressive PTI			3.8		100480	Type 49	
	Implant Progressive PTI			4.2		100480	Type 49	
	Implant Progressive PTI			5.0		100480	Type 49	
BIODENTA	Tissue level		NP			100402	Type 21	
	Tissue level		RP			100478	Type 47	
	Tissue level		WP			100478	Type 47	
BIOHORIZONS	External			3.5		100436	Type 30	
	External			4.0		100353	Type 01	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
BIOHORIZONS	External				5.0		100353	Type 01
	External				6.0		100353	Type 01
	Internal		Yellow				100440	Type 32
	Internal		Green				100431	Type 27
	Internal		Blue				100442	Type 33
	Internal Laser Lok				3.0		100425	Type 26
	Single Stage		Yellow				100440	Type 32
	Single Stage		Green				100431	Type 27
	Single Stage		Blue				100442	Type 33
	Tapered Internal		Yellow				100440	Type 32
	Tapered Internal		Green				100431	Type 27
	Tapered Internal		Blue				100442	Type 33
	BIOINFINITY	Junior Platform				3.2		100455
Standard Platform					4.8		100425	Type 26
BIOMET 3i	OSSEOTITE Certain		XP		4/5		100386	Type 15
	OSSEOTITE Certain		XP		5/6		100386	Type 15
	OSSEOTITE Tapered Certain		Purple				100472	Type 45
	OSSEOTITE Tapered Certain		Blue				100386	Type 15
	NanoTite Certain Microminiplant				3.25		100472	Type 45
	NanoTite External Hex		Purple				100353	Type 01
	NanoTite External Hex		Blue				100353	Type 01
	NanoTite External Hex		Yellow				100353	Type 01
	NanoTite External Hex		Green				100353	Type 01
	NanoTite External Hex				XP 3/4		100353	Type 01
	NanoTite External Hex				XP 4/5		100353	Type 01
	NanoTite External Hex				XP 5/6		100353	Type 01
	NanoTite External Hex MicroMiniplant		Purple				100353	Type 01
	NanoTite Tapered Certain		Purple				100472	Type 45
	NanoTite Tapered Certain		Blue				100386	Type 15
	NanoTite Tapered Certain		Yellow				100386	Type 15
	NanoTite Tapered Certain		Green				100386	Type 15
	NanoTite Tapered Certain Prevail		Purple				100472	Type 45
	NanoTite Tapered Certain Prevail		Blue				100386	Type 15

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
BIOMET 3i	NanoTite Tapered Certain Prevail		Yellow				100386	Type 15
	OSSEOTITE Certain MicroMiniplant			3.25			100472	Type 45
	OSSEOTITE External Hex		Purple				100353	Type 01
	OSSEOTITE External Hex		Blue				100353	Type 01
	OSSEOTITE External Hex		Yellow				100353	Type 01
	OSSEOTITE External Hex		Green				100353	Type 01
	OSSEOTITE External Hex			XP 3/4			100353	Type 01
	OSSEOTITE External Hex			XP 4/5			100353	Type 01
	OSSEOTITE External Hex			XP 5/6			100353	Type 01
	OSSEOTITE External Hex MicroMiniplant		Purple				100353	Type 01
	OSSEOTITE Tapered Certain		Yellow				100386	Type 15
	OSSEOTITE Tapered Certain		Green				100386	Type 15
	OSSEOTITE Tapered Certain Prevail		Purple				100472	Type 45
	OSSEOTITE Tapered Certain Prevail		Blue				100386	Type 15
	OSSEOTITE Tapered Certain Prevail		Yellow				100386	Type 15
	T3 Non-Platform Switched		Purple				100472	Type 45
	T3 Non-Platform Switched		Blue				100386	Type 15
	T3 Non-Platform Switched		Yellow				100386	Type 15
	T3 Non-Platform Switched		Green				100386	Type 15
	T3 Platform Switched		Purple				100472	Type 45
	T3 Platform Switched		Blue				100386	Type 15
	T3 Platform Switched		Yellow				100386	Type 15
	BIOTEC SRL	EN			All			100436
ER				All			100353	Type 01
IM				All			100431	Type 27
IR				All			100440	Type 32
KR				3.0			100455	Type 38
KW				3.6			100440	Type 32
SR				All			100378	Type 06
TWO				All			100440	Type 32
BIOTECH DENTAL	Kontakt K30			3.0			100581	Type 75 New
	Kontakt			3.6			100582	Type 76 New
	Kontakt			4.2			100582	Type 76 New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
BIOTECH DENTAL	Kontakt			4.8		100582	Type 76	New
	Kontakt			5.4		100582	Type 76	New
	Weego.FR			3.6		100583	Type 78	New
	Weego.FR			4.2		100583	Type 78	New
	Weego.FR			4.8		100583	Type 78	New
BONE SYSTEM	2P			5.0		100447	Type 36	
BRAINBASE			C			100478	Type 47	
	Mytis Arrow Implants		B			100350	Type 04	
	Mytis Arrow Implants		BW			100350	Type 04	
	Mytis Arrow Implants		E			100350	Type 04	
	Mytis Arrow Implants		EW			100350	Type 04	
BREIDENT	blueSKY					100480	Type 49	
BTI BIOTECHNOLOGY INSTITUTE	Externa		Universal			100353	Type 01	
	Externa		Universal plus			100353	Type 01	
	Interna		Universal			100431	Type 27	
	Interna		Universal plus			100431	Type 27	
	Interna		Wide			100431	Type 27	
	Tiny			2.5		100368	Type 09	
	Tiny			3.75		100368	Type 09	
BTLOCK	CV1			3.75		100480	Type 49	
	CV1			4.5		100480	Type 49	
	CV1			5.5		100480	Type 49	
	CV1			6.5		100480	Type 49	
	CV2			3.75		100480	Type 49	
	CV2			4.5		100480	Type 49	
	CV2			5.5		100480	Type 49	
	CV2			6.5		100480	Type 49	
	CV3			3.75		100480	Type 49	
	CV3			4.5		100480	Type 49	
	CV3			5.5		100480	Type 49	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
BTLOCK	CV3			6.5		100480	Type 49	
CAMLOG	Camlog Cylinder-line			3.3		100432	Type 28	
	Camlog Cylinder-line			3.8		100411	Type 23	
	Camlog Cylinder-line			4.3		100411	Type 23	
	Camlog Cylinder-line			5.0		100412	Type 24	
	Camlog Cylinder-line			6.0		100412	Type 24	
	Camlog Root-line Promote			3.8		100411	Type 23	
	Camlog Root-line Promote			4.3		100411	Type 23	
	Camlog Root-line Promote			5.0		100412	Type 24	
	Camlog Root-line Promote			6.0		100412	Type 24	
	Camlog Screw-line			3.3		100432	Type 28	
	Camlog Screw-line			3.8		100411	Type 23	
	Camlog Screw-line			4.3		100411	Type 23	
	Camlog Screw-line			5.0		100412	Type 24	
	Camlog Screw-line			6.0		100412	Type 24	
	Conelog Screw-line			3.3		100557	Type 65	
	Conelog Screw-line			3.8		100541	Type 58	
	Conelog Screw-line			4.3		100541	Type 58	
	Conelog Screw-line			5.0		100380	Type 07	
	iSy (direct to implant)			All		100497	Type 57	
	iSy (to implant base)			All		100446	Type 35	
CHANG GUNG	ComMed			3.5		100481	Type 50	
	ComMed			4.0		100481	Type 50	
	ComMed			4.5		100380	Type 07	
	ComMed			5.0		100380	Type 07	
	ComMed			5.5		100380	Type 07	
	ComMed			6.0		100380	Type 07	
CHAORUM	External			3.3		100353	Type 01	
	Internal			4.0		100478	Type 47	
	Internal			4.8		100478	Type 47	
COCHLEAR	Baha & Vistafix Flange Fixture Hexagonal Connection			3.0		100368	Type 09	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
COCHLEAR	Baha & Vistafix Flange Fixture Hexagonal Connection			4.0		100368	Type 09	
	Baha BI 300			3.0		100436	Type 30	
	Baha BI 300			4.0		100436	Type 30	
	Baha BIA 300			6.0		100487	Type 55	
	Baha BIA 300			9.0		100487	Type 55	
	Baha BIA 400			6.0		100487	Type 55	
	Baha BIA 400			8.0		100487	Type 55	
	Baha BIA 400			10.0		100487	Type 55	
	Baha BIA 400			12.0		100487	Type 55	
	Baha BIA 400			14.0		100487	Type 55	
	Vistafix VXI 300			3.0		100436	Type 30	
	Vistafix VXI 300			4.0		100436	Type 30	
COMPRESS				4.0		100355	Type 02	
CORTEX	Classix			3.3		100425	Type 26	
	Classix			3.8		100425	Type 26	
	Classix			4.2		100425	Type 26	
	Classix			5.0		100425	Type 26	
	Classix			6.0		100425	Type 26	
	Dynamix			3.3		100425	Type 26	
	Dynamix			3.8		100425	Type 26	
	Dynamix			4.2		100425	Type 26	
	Dynamix			5.0		100425	Type 26	
	Dynamix			6.0		100425	Type 26	
	Saturn			3.3		100425	Type 26	
	Saturn			3.8		100425	Type 26	
	Saturn			4.2		100425	Type 26	
	Saturn			5.0		100425	Type 26	
	Saturn			6.0		100425	Type 26	
COWELLMEDI	Atlas External			3.3		100355	Type 02	
	Atlas External			4.1		100353	Type 01	
	Atlas External			5.0		100353	Type 01	
	Atlas Submerged			4.0		100402	Type 21	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
C-TECH	EL				3.5		100425	Type 26
	EL				4.3		100425	Type 26
	EL				5.1		100425	Type 26
DENTAL TECH			WP				100478	Type 47
DENTAURUM	Tiologic		S				100472	Type 45
	Tiologic		M				100472	Type 45
	Tiologic		L				100472	Type 45
DENTEGRIS	Straight & Tapered				3.3		100425	Type 26
	Straight & Tapered				3.75		100440	Type 32
	Straight & Tapered				4.5		100440	Type 32
	Straight & Tapered				4.75		100440	Type 32
	Straight & Tapered				5.5		100478	Type 47
DENTIS DENTAL	Internal				3.7		100350	Type 04
	Internal				4.3		100350	Type 04
	Internal				4.8		100350	Type 04
	OneQ S-clean Narrow (SLA surface/Submerged)				3.0		100404	Type 22
	OneQ S-clean Narrow (SLA surface/Submerged)				3.3		100404	Type 22
	OneQ S-clean Regular (SLA surface/Submerged)				3.7		100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				4.2		100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				4.7		100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				5.2		100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				6.0		100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				7.0		100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				8.0		100372	Type 11
	Save Submerged				5.5		100380	Type 07
	Save Submerged				6.0		100380	Type 07
	S-clean Mini (Submerged)				3.7		100372	Type 11
	S-clean Regular (Submerged)				4.1		100372	Type 11
	S-clean Regular (Submerged)				4.3		100372	Type 11
	S-clean Wide (Submerged)				4.8		100372	Type 11

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
DENTIUM	Genoss			3.0			100554	Type 63
	Genoss			3.6			100402	Type 21
	Genoss			4.3			100402	Type 21
	Genoss			5.0			100402	Type 21
	Implantium		Yellow 3.6	3.4			100378	Type 06
	Implantium		Green 4.0	3.8			100378	Type 06
	Implantium		Blue 4.5	4.3			100378	Type 06
	Simple Line			3.8			100380	Type 07
	Simple Line			4.3			100380	Type 07
	Simple Line			4.8			100380	Type 07
	Super Line			4.0			100380	Type 07
	Super Line			7.0			100380	Type 07
DENTSPLY IMPLANTS	ANKYLOS C/X Implants		A				100388	Type 16
	ANKYLOS C/X Implants		B				100388	Type 16
	ANKYLOS C/X Implants		C				100388	Type 16
	ASTRA TECH Implant system MicroThread			3.5			100376	Type 05
	ASTRA TECH Implant system MicroThread			4.0			100378	Type 06
	ASTRA TECH Implant system MicroThread			4.5			100380	Type 07
	ASTRA TECH Implant system MicroThread			5.0			100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Aqua			3.5 S			100455	Type 38
	ASTRA TECH Implant system OsseoSpeed Aqua			4.0 S			100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25212	Green	3.0 S	8		100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25213	Green	3.0 S	9		100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25214	Green	3.0 S	11		100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25215	Green	3.0 S	13		100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25216	Green	3.0 S	15		100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25221	Purple	3.6 S	6		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25222	Purple	3.6 S	8		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25223	Purple	3.6 S	9		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25224	Purple	3.6 S	11		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25225	Purple	3.6 S	13		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25226	Purple	3.6 S	15		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25227	Purple	3.6 S	17		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25231	Yellow	4.2 S	6		100480	Type 49

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
DENTSPLY IMPLANTS	ASTRA TECH Implant system OsseoSpeed EV	25232	Yellow		4.2 S 8		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25233	Yellow		4.2 S 9		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25234	Yellow		4.2 S 11		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25235	Yellow		4.2 S 13		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25236	Yellow		4.2 S 15		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25237	Yellow		4.2 S 17		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25241	Blue		4.8 S 6		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25242	Blue		4.8 S 8		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25243	Blue		4.8 S 9		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25244	Blue		4.8 S 11		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25245	Blue		4.8 S 13		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25246	Blue		4.8 S 15		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25247	Blue		4.8 S 17		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25251	Brown		5.4 S 6		100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25252	Brown		5.4 S 8		100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25253	Brown		5.4 S 9		100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25254	Brown		5.4 S 11		100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25255	Brown		5.4 S 13		100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25256	Brown		5.4 S 15		100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25262	Yellow		4.2 C 8		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25263	Yellow		4.2 C 9		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25264	Yellow		4.2 C 11		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25265	Yellow		4.2 C 13		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25266	Yellow		4.2 C 15		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25267	Yellow		4.2 C 17		100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25272	Blue		4.8 C 8		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25273	Blue		4.8 C 9		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25274	Blue		4.8 C 11		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25275	Blue		4.8 C 13		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25276	Blue		4.8 C 15		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25277	Blue		4.8 C 17		100478	Type 47
	ASTRA TECH Implant system OsseoSpeed Lilac				4.5		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Lilac				5.0		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Lilac				5.0S		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Profile EV	25452	Yellow		4.2 PC 8		100370	Type 10

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
DENTSPLY IMPLANTS	ASTRA TECH Implant system OsseoSpeed Profile EV	25453	Yellow	4.2 PC	9		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25454	Yellow	4.2 PC	11		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25455	Yellow	4.2 PC	13		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25456	Yellow	4.2 PC	15		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25457	Yellow	4.2 PC	17		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25458	Blue	4.8 PC	8		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25459	Blue	4.8 PC	9		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25460	Blue	4.8 PC	11		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25461	Blue	4.8 PC	13		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25462	Blue	4.8 PC	15		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25463	Blue	4.8 PC	17		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25440	Yellow	4.2 PS	8		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25441	Yellow	4.2 PS	9		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25442	Yellow	4.2 PS	11		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25443	Yellow	4.2 PS	13		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25444	Yellow	4.2 PS	15		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25445	Yellow	4.2 PS	17		100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25446	Blue	4.8 PS	8		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25447	Blue	4.8 PS	9		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25448	Blue	4.8 PS	11		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25449	Blue	4.8 PS	13		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25450	Blue	4.8 PS	15		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25451	Blue	4.8 PS	17		100440	Type 32
	ASTRA TECH Implant system OsseoSpeed TX Aqua			3.5S			100455	Type 38
	ASTRA TECH Implant system OsseoSpeed TX Aqua			4.0S			100455	Type 38
	ASTRA TECH Implant system OsseoSpeed TX Lilac			4.5			100380	Type 07
	ASTRA TECH Implant system OsseoSpeed TX Lilac			5.0			100380	Type 07
	ASTRA TECH Implant system OsseoSpeed TX Lilac			5.0 S			100380	Type 07
	ASTRA TECH Implant system OsseoSpeed TX Profile			4.5			100378	Type 06
	ASTRA TECH Implant system OsseoSpeed TX Profile			5.0			100378	Type 06
	ASTRA TECH Implant system OsseoSpeed TX Profile			5.0 S			100378	Type 06
	ASTRA TECH Implant system OsseoSpeed TX Yellow			3.0 S			100481	Type 50
	ASTRA TECH Implant system OsseoSpeed Yellow			3.0 S			100481	Type 50
	XiVe		TG				100366	Type 08
	XiVe			3.0			100404	Type 22

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
DENTSPLY IMPLANTS	XiVe			3.4		100402	Type 21	
	XiVe			3.8		100472	Type 45	
	XiVe			4.5		100472	Type 45	
	XiVe			5.5		100472	Type 45	
DIO	External FSN			3.3		100355	Type 02	
	External FSN			3.75		100353	Type 01	
	External FTN			5.0		100357	Type 03	
	SM Submerged			3.8		100480	Type 49	
	SM Submerged			4.5		100478	Type 47	
	SM Submerged			5.3		100478	Type 47	
DTI				3.5		100497	Type 57	
				4.2		100380	Type 07	
DYNA DENTAL	Dyna Helix DC			3.2		100431	Type 27	
	Dyna Helix DC			3.6		100431	Type 27	
	Dyna Helix DC			4.2		100431	Type 27	
	Dyna Helix DC			5.0		100431	Type 27	
	Dyna Helix ST			3.6		100431	Type 27	
	Dyna Helix ST			4.2		100431	Type 27	
	Dyna Helix ST			5.0		100431	Type 27	
ECKERMANN	All Spiral					100390	Type 17	
	Hexagono Ext Hex					100353	Type 01	
	Triplo					100480	Type 49	
EMFILS	Colosso IH ECI Evolution			3.3		100372	Type 11	
	Colosso IH ECI Evolution			4.0		100372	Type 11	
	Colosso IH ECI Evolution			5.0		100372	Type 11	
	ICM-4010			4.0		100484	Type 53	
	PAT5-4010			4.0		100484	Type 53	
ENTIFIC	Flange Fixture					100353	Type 01	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
EUROTEKNIKA	Aesthetica+2		3.6		4.2		100544	Type 61
	Aesthetica+2		3.6		4.8		100544	Type 61
	Aesthetica+2		4.1		4.2		100544	Type 61
	Aesthetica+2		4.1		4.8		100544	Type 61
	Aesthetica+2		4.8		4.8		100544	Type 61
	Aesthetica+2		4.8		6.5		100544	Type 61
	Naturall				3.5		100455	Type 38
	Naturall				4.0		100378	Type 06
	Naturall				4.5		100378	Type 06
	Naturall				5.0		100378	Type 06
	Uneva+				3.6		100353	Type 01
	Uneva+				4.1		100353	Type 01
	Uneva+				4.8		100353	Type 01
	Uneva+				6.0		100353	Type 01
FAIR IMPLANT	FairTwo				3.5		100484	Type 53
	FairTwo				4.2		100484	Type 53
	FairTwo				5.0		100378	Type 06
	FairTwo				6.0		100378	Type 06
GC IMPLANT	Aadva		Narrow				100455	Type 38
	Aadva		Regular				100541	Type 58
	Aadva		Wide				100541	Type 58
	External Plus				3.0		100402	Type 21
	External Plus				3.4		100402	Type 21
	External Plus				3.8		100353	Type 01
	External Plus				5.0		100353	Type 01
	GENESIO						100402	Type 21
	Internal Plus				3.0		100402	Type 21
	Internal Plus				3.4		100402	Type 21
	SETIO						100353	Type 01
GHIMAS	PFON EVO SERIES				3.6		100455	Type 38
	PFON EVO SERIES				4.0		100455	Type 38
	PFON EVO SERIES				5.0		100380	Type 07

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
HIOSSEN	ETIII		Mini		3.5		100376	Type 05
	ETIII				4.0		100478	Type 47
	ETIII				4.5		100478	Type 47
	ETIII				5.0		100478	Type 47
	ETIII		UW		6.0		100378	Type 06
	ETIII		UW		7.0		100378	Type 06
	HGII		Mini				100376	Type 05
	HGII		Regular				100378	Type 06
	HGII		Wide				100378	Type 06
	HGII		Ultra Wide				100378	Type 06
	HGIII		Mini				100376	Type 05
	HSII		Regular				100350	Type 04
	HSII		Wide				100350	Type 04
IBS INNO BIO SURG	Magicore				5.5		100455	Type 38
IDI EVOLUTION	Vite				3.75		100398	Type 18
	Vite				4.0		100398	Type 18
IMMEDIATELOAD	IM40 Immediate Conico				4.0		100480	Type 49
IMPLANT DIRECT	Legacy1		Blue				100455	Type 38 New
	Legacy1		Green				100440	Type 32 New
	Legacy1		Purple				100431	Type 27 New
	Legacy1		Yellow				100442	Type 33 New
	Legacy2		Blue				100455	Type 38 New
	Legacy2		Green				100440	Type 32 New
	Legacy2		Purple				100431	Type 27 New
	Legacy2		Yellow				100442	Type 33 New
	Legacy3		Blue				100455	Type 38
	Legacy3		Green				100440	Type 32 New
	Legacy3		Purple				100431	Type 27 New
	Legacy3		Yellow				100442	Type 33 New
	Legacy4		Blue				100455	Type 38 New
	Legacy4		Green				100440	Type 32 New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New	
IMPLANT DIRECT	Legacy4		Purple				100431	Type 27	New
	Legacy4		Yellow				100442	Type 33	New
	ScrewPlant			3.7			100425	Type 26	
	ScrewPlant			4.7			100431	Type 27	
	ScrewPlant			5.7			100431	Type 27	
IMTEC	Hexed-head			3.75			100353	Type 01	
	Hexed-head			4.75			100353	Type 01	
INTRA-LOCK SYSTEM	Internal		SQ				100378	Type 06	
	Internal		Standard				100478	Type 47	
	Internal		Wide				100478	Type 47	
IRES	IMAX NHSE			3.3			100355	Type 02	
	IMAX NHSE			3.7			100353	Type 01	
	IMAX NHSE			4.1			100353	Type 01	
	IMAX NHSE			5.2			100353	Type 01	
	IMAX NHSI			3.7			100431	Type 27	
	IMAX NHSI			4.1			100431	Type 27	
	IMAX NHSI			5.2			100431	Type 27	
	Shapeone S1B			3.1			100431	Type 27	
	Shapeone S1B			4.1			100431	Type 27	
	Shapeone S1B			4.7			100431	Type 27	
	Shapeone S1EH			3.1			100353	Type 01	
	Shapeone S1EH			4.1			100353	Type 01	
	Shapeone S1EH			4.7			100353	Type 01	
	Shapeone S1T			3.7			100350	Type 04	
	Shapeone S1T			4.1			100350	Type 04	
	Shapeone S1T			4.7			100350	Type 04	
	Shapeone S1TN			3.7			100431	Type 27	
	Shapeone S1TN			4.1			100431	Type 27	
	Shapeone S1TN			4.7			100431	Type 27	
	Volution SVB			3.7			100431	Type 27	
	Volution SVB			4.1			100431	Type 27	
Volution SVB			4.7			100431	Type 27		

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
ITERUM	Conical			3.75		100425	Type 26	
JDENTALCARE	JDEvolution			3.2		100455	Type 38	
	JDEvolution			3.7		100440	Type 32	
	JDEvolution			4.3		100440	Type 32	
	JDEvolution			5.0		100440	Type 32	
JMM	POI 2 piece			3.2		100398	Type 18	
	POI 2 piece			3.7		100398	Type 18	
	POI 2 piece			4.2		100398	Type 18	
	POIEX			3.7		100386	Type 15	
	POIEX			4.2		100480	Type 49	
	POIEX			4.7		100478	Type 47	
	POIEX			5.2		100478	Type 47	
KENTEC	SB1			3.3		100378	Type 06	
	SB2			4.2		100378	Type 06	
KEYSTONE DENTAL	Genesis			3.8		100431	Type 27	
	Genesis			4.5		100442	Type 33	
	Genesis			5.5		100442	Type 33	
	Genesis			6.5		100442	Type 33	
	Prima Connex		SD	3.5		100431	Type 27	
	Prima Connex		RD	4.1		100442	Type 33	
	Prima Connex		WD	5.0		100442	Type 33	
	Renova RDI			4.5		100431	Type 27	
	Renova SDI			3.75		100425	Type 26	
	Restore			3.3		100436	Type 30	
	Restore			3.75		100353	Type 01	
	Restore			4.0		100353	Type 01	
	Restore			5.0		100357	Type 03	
	Restore			6.0		100357	Type 03	
KLOCKNER IMPLANT SYSTEM	MV VEGA Implant			3.0		100497	Type 57	
	NV VEGA Implant			3.5		100497	Type 57	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
KLOCKNER IMPLANT SYSTEM	RV VEGA Implant				4.0		100425	Type 26
	RV VEGA Implant				4.5		100425	Type 26
	Essential Cone				3.5		100480	Type 49
	Essential Cone				4.0		100480	Type 49
	Essential Cone				4.5		100480	Type 49
	Essential Cone				4.8		100480	Type 49
	KL Implant NP				3.3		100355	Type 02
	KL Implant NP				3.7		100353	Type 01
	KL Implant NP				4.2		100353	Type 01
	KL Implant WP				4.7		100357	Type 03
	S3M						100436	Type 30
SK2						100436	Type 30	
KYOCERA	FINESIA BL		NP				100545	Type 62
	FINESIA BL		RP				100431	Type 27
	FINESIA BL		WP				100350	Type 04
	FINESIA TL		RP				100350	Type 04
	FINESIA TL		WP				100350	Type 04
LASAK	BioniQ, QN (Q-Lock Narrow)						100404	Type 22
	BioniQ, QR (Q-Lock Regular)						100497	Type 57
	Impladent				2.9		100355	Type 02
	Impladent				3.7		100447	Type 36
	Impladent				5.0		100448	Type 37
LEADER	Implus EH				3.3		100353	Type 01
	Implus EH				3.75		100353	Type 01
	Implus EH				5.0		100353	Type 01
	Implus IH				3.3		100425	Type 26
	Tixos IH				3.75		100425	Type 26
	Tixos IH				4.5		100425	Type 26
	Tixos IH				5.5		100425	Type 26
	Tixos ITMC				3.75		100497	Type 57
	Tixos ITMC				4.5		100497	Type 57
	Tixos ITMC				5.5		100497	Type 57

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
MEDENTIKA	M-Implant®			3.5		100455	Type 38	
	M-Implant®			4.0		100455	Type 38	
	M-Implant®			4.5		100455	Type 38	
	M-Implant®			5.0		100455	Type 38	
MEDICAL INSTINCT	Bone Trust			3.4		100466	Type 44	
	Bone Trust			4.0		100466	Type 44	
	Bone Trust			5.0		100466	Type 44	
	Bone Trust Plus Cone			3.4		100466	Type 44	
	Bone Trust Plus Cone			4.0		100466	Type 44	
MEDIGMA TECHNOLOGIES	Fix 1					100390	Type 17	
	Fix-A-Dent					100480	Type 49	
	FixTite					100480	Type 49	
	FixTite-S					100480	Type 49	
MEGAGEN	AnyOne			3.5		100380	Type 07	
	AnyOne			4.0		100380	Type 07	
	AnyOne External			4.5		100350	Type 04	
	AnyRidge					100431	Type 27	
	EZ Plus Internal		Small			100455	Type 38	
	EZ Plus Internal		Regular			100380	Type 07	
	EZ Plus Internal		Wide			100380	Type 07	
	EZ Plus Internal			4.0		100380	Type 07	
	EZ Plus Internal			4.5		100380	Type 07	
	EZ Plus Internal			5.0		100380	Type 07	
	Mini Internal			3.0		100481	Type 50	
	Resque External			7.0		100357	Type 03	
MICRODENT	Duraplant 2.2			4.1		100370	Type 10	
	Microdent system					100353	Type 01	
	MK			5.6		100357	Type 03	
	MS			3.5		100436	Type 30	
	MST			3.5		100436	Type 30	
	Trylogic			3.8		100431	Type 27	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
MICRODENT	Trylogic			4.3			100431	Type 27
	Universal			4.1			100353	Type 01
	Universal			5.1			100353	Type 01
MIS	C1 Conical Connection		NP				100455	Type 38
	C1 Conical Connection		SP				100480	Type 49
	C1 Conical Connection		WP				100584	Type 79
	V3 Conical Connection		NP				100455	Type 38
	V3 Conical Connection		SP				100480	Type 49
	Internal Hexagon (seven, biocom)				3.75		100440	Type 32
	Internal Octagon (mistral)				4.75		100350	Type 04
Seven		Narrow		3.3		100554	Type 63	
MOZO-GRAU	MG-Inhex		Standard				100455	Type 38
	MG-Inhex		Maxi				100380	Type 07
	MG-Osseous		Mini				100353	Type 01
	MG-Osseous		Standard				100353	Type 01
	MG-Osseous		Maxi				100357	Type 03
ND IMPLANTS	Niko			3.6			100431	Type 27
	Niko			4.6			100380	Type 07
NEMRIS	Aesthura Classic			3.75			100446	Type 35
	Aesthura Classic			4.75			100446	Type 35
NEOBIOTECH	IS Fixture (RBM Surface)			4.0			100376	Type 05
	IS Fixture (RBM Surface)			4.5			100376	Type 05
	IS Fixture (RBM Surface)			4.8			100376	Type 05
	IS S-Wide Fixture (RBM Surface)			5.5			100376	Type 05
	IS S-Wide Fixture (RBM Surface)			6.0			100376	Type 05
	IS S-Wide Fixture (RBM Surface)			7.0			100376	Type 05
	IS-II active Fixture (S.L.A. Surface)			3.5			100376	Type 05
	IS-II active Fixture (S.L.A. Surface)			4.0			100376	Type 05
	IS-II active Fixture (S.L.A. Surface)			4.5			100376	Type 05
	IS-II active Fixture (S.L.A. Surface)			5.0			100376	Type 05

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
NEOBIOTECH	IS-II active S-Wide Fixture (S.L.A. Surface)				5.5		100376	Type 05
	IS-II active S-Wide Fixture (S.L.A. Surface)				6.0		100376	Type 05
	IS-II active S-Wide Fixture (S.L.A. Surface)				7.0		100376	Type 05
	IS-II Fixture (RBM Surface)				3.5		100376	Type 05
	IS-II Fixture (RBM Surface)				4.0		100376	Type 05
	IS-II Fixture (RBM Surface)				4.5		100376	Type 05
	IS-II Fixture (RBM Surface)				5.0		100376	Type 05
	IS-II S-Wide Fixture (RBM Surface)				5.5		100376	Type 05
	IS-II S-Wide Fixture (RBM Surface)				6.0		100376	Type 05
	IS-II S-Wide Fixture (RBM Surface)				7.0		100376	Type 05
	IS-III active Fixture				3.5		100376	Type 05
	IS-III active Fixture				4.0		100376	Type 05
	IS-III active Fixture				4.5		100376	Type 05
	IS-III active Fixture				5.0		100376	Type 05
	IS-III active S-Wide Fixture				5.5		100376	Type 05
	IS-III active S-Wide Fixture				6.0		100376	Type 05
	IS-III active S-Wide Fixture				7.0		100376	Type 05
	Neoplant External				4.0		100353	Type 01
	Neoplant External				5.0		100353	Type 01
NEODENT	Alvim CM				3.5		100388	Type 16
	Alvim CM				4.3		100388	Type 16
	Alvim CM				5.0		100388	Type 16
	Alvim II Plus				3.75		100402	Type 21
	Alvim II Plus				4.3		100402	Type 21
	Alvim II Plus				5.0		100402	Type 21
	Alvim Ti				4.3		100353	Type 01
	Drive CM				4.3		100479	Type 48
	Drive CM				5.0		100479	Type 48
	Titamax CM				3.5		100388	Type 16
	Titamax CM				3.75		100388	Type 16
	Titamax CM				4.0		100388	Type 16
	Titamax CM				5.0		100388	Type 16
	Titamax GT				3.7		100368	Type 09
	Titamax GT				4.0		100368	Type 09

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
NEODENT	Titamax GT				4.5		100368	Type 09
	Titamax II Plus				3.75		100402	Type 21
	Titamax II Plus				5.0		100402	Type 21
	Titamax Ti				3.75		100353	Type 01
	Titamax Ti				4.0		100353	Type 01
	Titamax Ti				4.5		100357	Type 03
	Titamax Ti				5.0		100357	Type 03
NEOSS	Neoss Implant				3.25		100455	Type 38
	ProActive Straight				3.5		100398	Type 18
	ProActive Straight				4.0		100398	Type 18
	ProActive Straight				4.5		100398	Type 18
	ProActive Straight				5.0		100398	Type 18
	ProActive Straight				5.5		100398	Type 18
	ProActive Straight				6.0		100398	Type 18
	ProActive Tapered				3.5		100398	Type 18
	ProActive Tapered				4.0		100398	Type 18
	ProActive Tapered				4.5		100398	Type 18
	ProActive Tapered				5.0		100398	Type 18
	ProActive Tapered				5.5		100398	Type 18
	ProActive Tapered				6.0		100398	Type 18
NIPPON PISTON RING	1 Stage Screw RH				3.3		100398	Type 18
	1 Stage Screw RH				4.0		100398	Type 18
	1 Stage Screw RL				3.3		100398	Type 18
	1 Stage Screw RL				4.0		100398	Type 18
	1 Stage Screw WH				5.0		100398	Type 18
	1 Stage Screw WL				5.0		100398	Type 18
	2 Stage Cylinder				3.3		100398	Type 18
	2 Stage Cylinder				4.0		100398	Type 18
	2 Stage Cylinder				5.0		100398	Type 18
	2 Stage Screw				3.3		100398	Type 18
	2 Stage Screw				4.0		100398	Type 18
	2 Stage Screw				5.0		100398	Type 18

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
NOBEL BIOCARE	Brånemark system		NP				100355	Type 02
	Brånemark system		RP				100353	Type 01
	Brånemark system		WP				100357	Type 03
	Brånemark system Novum						100353	Type 01
	Brånemark system Shorty		RP				100353	Type 01
	Brånemark system Shorty		WP				100357	Type 03
	NobelActive		NP	3.5			100543	Type 60
	NobelActive		RP	4.3			100544	Type 61
	NobelActive		RP	5.0			100544	Type 61
	NobelActive		WP	5.5			100350	Type 04
	NobelActive			3.0			100545	Type 62
	NobelDirect		NP				100366	Type 08
	NobelDirect		RP				100368	Type 09
	NobelDirect		WP				100368	Type 09
	NobelParallel Conical Connection		NP	3.75			100543	Type 60
	NobelParallel Conical Connection		RP	4.3			100544	Type 61
	NobelParallel Conical Connection		RP	5.0			100544	Type 61
	NobelParallel Conical Connection		WP	5.5			100350	Type 04
	NobelReplace Conical Connection		NP	3.5			100543	Type 60
	NobelReplace Conical Connection		RP	4.3			100544	Type 61
	NobelReplace Conical Connection		RP	5.0			100544	Type 61
	NobelReplace Conical Connection PMC		NP				100543	Type 60
	NobelReplace Conical Connection PMC		RP				100544	Type 61
	NobelReplace Platform Shift		NP				100359	Type 12
	NobelReplace Platform Shift		RP				100361	Type 13
	NobelReplace Platform Shift		WP				100361	Type 13
	NobelReplace Straight Groovy		NP				100359	Type 12
	NobelReplace Straight Groovy		RP				100361	Type 13
	NobelReplace Straight Groovy		WP				100361	Type 13
	NobelReplace Tapered Groovy		NP				100359	Type 12
	NobelReplace Tapered Groovy		RP				100361	Type 13
	NobelReplace Tapered Groovy		WP				100361	Type 13
	NobelReplace Tapered Groovy			6.0			100462	Type 40
	NobelSpeedy Groovy		NP				100355	Type 02
	NobelSpeedy Groovy		RP				100353	Type 01

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
NOBEL BIOCARE	NobelSpeedy Groovy		WP				100357	Type 03
	NobelSpeedy Replace		NP				100359	Type 12
	NobelSpeedy Replace		RP				100361	Type 13
	NobelSpeedy Replace		WP				100361	Type 13
	NobelSpeedy Replace			6.0			100462	Type 40
	NobelSpeedy Shorty		RP				100353	Type 01
	NobelSpeedy Shorty		WP				100357	Type 03
	Replace Select Straight		NP				100359	Type 12
	Replace Select Straight		RP				100361	Type 13
	Replace Select Straight		WP				100361	Type 13
	Replace Select Tapered		NP				100359	Type 12
	Replace Select Tapered		RP				100361	Type 13
	Replace Select Tapered		WP				100361	Type 13
	Replace Select Tapered			6.0			100462	Type 40
	Replace Select Tapered PMC		NP				100359	Type 12
	Replace Select Tapered PMC		RP				100361	Type 13
	Replace Select Tapered PMC		WP				100361	Type 13
	Replace Select Tapered PMC			6.0			100462	Type 40
	Replace Select TC		NP				100359	Type 12
	Replace Select TC		RP				100361	Type 13
NUCLEOSS	Nucleoss T3			3.4			100376	Type 05
	Nucleoss T3			3.8			100376	Type 05
	Nucleoss T3			4.2			100376	Type 05
	Nucleoss T3			5.0			100376	Type 05
	Nucleoss T4			3.4			100376	Type 05
	Nucleoss T4			3.8			100376	Type 05
	Nucleoss T4			4.2			100376	Type 05
	Nucleoss T5			5.0			100380	Type 07
OCO BIOMEDICAL	Engage		3.0	3.2			100455	Type 38
	Engage		3.5	4.0			100425	Type 26
	Engage		4.5	5.0			100425	Type 26
	Engage		4.5	6.0			100431	Type 27
	ERI			4.0			100380	Type 07

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
OCO BIOMEDICAL	ERI				5.0		100380	Type 07
	Macro				6.0		100478	Type 47
	TSI				3.25		100446	Type 35
	TSI				4.0		100380	Type 07
	TSI				5.0		100380	Type 07
OSSEOLINK					3.5		100402	Type 21
					4.0		100402	Type 21
					5.0		100402	Type 21
OSSTEM	GSII		Mini				100376	Type 05
	GSII		Regular				100378	Type 06
	GSII		Wide				100378	Type 06
	GSII		Ultra Wide				100378	Type 06
	GSIII		Mini				100376	Type 05
	GSIII		Regular				100378	Type 06
	SSI		Regular				100350	Type 04
	SSI		Wide				100350	Type 04
	SSII		Regular				100350	Type 04
	SSII		Wide				100350	Type 04
	SSIII		Regular				100350	Type 04
	SSIII		Wide				100350	Type 04
	TSIII CA				3.5		100376	Type 05
	TSIII SA				3.0		100372	Type 11
	TSIII SA				4.0		100478	Type 47
	TSIII SA				4.5		100478	Type 47
	TSIII SA				5.0		100478	Type 47
	TSIII SA Fixture – NoMount				4.5		100478	Type 47
	TSIII SA Ultra-Wide				6.0		100350	Type 04
	TSIII SA Ultra-Wide				7.0		100350	Type 04
	USII		Mini				100355	Type 02
	USII		Regular				100353	Type 01
	USII		Wide				100357	Type 03
	USII		Ultra Wide				100357	Type 03
	USIII		Mini				100355	Type 02

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
OSSTEM	USIII		Regular				100353	Type 01
	USIII		Wide				100357	Type 03
	USIII		Ultra Wide				100357	Type 03
	USIV Conical		Mini				100366	Type 08
OSTEOPLANT	Two-piece						100370	Type 10
OSTEO-TI	Combi Dental Implant				3.25		100425	Type 26
	Combi Dental Implant				3.75		100425	Type 26
	Combi Dental Implant				4.5		100425	Type 26
	Combi Dental Implant				5.5		100425	Type 26
OTICON MEDICAL	Ponto Implant				3.0		100368	Type 09
	Ponto Implant				4.0		100368	Type 09
	Ponto Implant		Wide				100368	Type 09
PALTOP	PALTOP		NP				100455	Type 38
	PALTOP		SP				100480	Type 49
PHIBO	Aurea EVO		NP				100554	Type 63
	Aurea EVO		RP				100425	Type 26
	Aurea EVO		WP				100425	Type 26
	Aurea		NP				100554	Type 63
	Aurea		RP				100431	Type 27
	Aurea		WP				100440	Type 32
	BNT				S3		100353	Type 01
	BNT				S4		100353	Type 01
	BNT				S5		100353	Type 01
	TSA				S3		100455	Type 38
	TSA				S4		100431	Type 27
	TSA				S5		100431	Type 27
	TSA Advance				S3		100455	Type 38
	TSA Advance				S4		100431	Type 27
	TSA Advance				S5		100431	Type 27
	TSH				S3		100353	Type 01

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
PHIBO	TSH			S4		100353	Type 01	
	TSH			S5		100353	Type 01	
P-I BRÄNEMARK PHILOSOPHY	Internal Grip Amplified			3.5		100554	Type 63	
	Internal Grip Amplified			4.3		100482	Type 51	
	Internal Grip Amplified			5.1		100482	Type 51	
	Internal Grip External Hex			3.5		100455	Type 38	
	Internal Grip External Hex			3.75		100398	Type 18	
	Internal Grip External Hex			4.1		100398	Type 18	
	Morse Taper			3.3		100402	Type 21	
	Morse Taper			3.75		100402	Type 21	
	Morse Taper			4.8		100402	Type 21	
PLAN 1 HEALTH	GTB			3.3		100497	Type 57	
	GTB			3.6		100497	Type 57	
	GTB			4.3		100497	Type 57	
	Hexcel EHC, External Hex					100353	Type 01	
	Hexcel IHC, Internal Hex					100425	Type 26	
	Hexcel IHT, Internal Hex					100425	Type 26	
PLATON	EL Implant Narrow			3.7		100543	Type 60	
	EL Implant Regular			5.1		100478	Type 47	
	Type I					100350	Type 04	
	Type II					100350	Type 04	
	Type IV Pro					100478	Type 47	
PRODENT ITALIA	Perform SM			3.6		100480	Type 49	
	Perform SM			4.0		100480	Type 49	
	Perform SM			4.5		100380	Type 07	
	Perform SM			5.0		100380	Type 07	
PW PLUS	Dif			3.8		100480	Type 49	New
	Dif			4.2		100380	Type 07	New
REVIOS				3.8		100425	Type 26	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
REVIOS				4.3		100425	Type 26	
				5.0		100425	Type 26	
REX IMPLANTS	Rex TL			1.8		100350	Type 04	
	Rex TL			2.9		100480	Type 49	
S.I.N IMPLANT SYSTEM	External Hex			3.4		100353	Type 01	
	External Hex			4.1		100353	Type 01	
	External Hex			5.0		100353	Type 01	
	Internal Hex			3.8		100386	Type 15	
	Internal Hex			4.5		100386	Type 15	
	Internal Hex			5.5		100386	Type 15	
	Strong SW HI			4.5		100479	Type 48	
	Strong SW Morse			4.5		100484	Type 53	
SDS	SDS 2.0			3.8		100589	Type 69	
SERVO-DENTAL	Trias			3.8		100378	Type 06	
SHINHUNG	Luna		NP	3.5		100376	Type 05	
	Luna		RP	4.0		100478	Type 47	
	Luna		RP	4.5		100478	Type 47	
	Luna		RP	5.0		100478	Type 47	
	Sola		NP	3.5		100355	Type 02	
	Sola		RP	4.0		100353	Type 01	
	Sola		RP	4.5		100353	Type 01	
	Sola		WP	5.0		100357	Type 03	
	Stella		RP	4.0		100544	Type 61	
	Stella		RP	4.5		100544	Type 61	
	Stella		WP	4.5		100544	Type 61	
	Stella		WP	5.0		100544	Type 61	
SHOFU	BFIMP30-OSTL			3.0		100545	Type 62	
	BFIMP35-OSTL			3.5		100404	Type 22	
	BFIMP40-OSTL			4.0		100543	Type 60	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
SIC IMPLANT SYSTEM	Ace			3.4		100466	Type 44	
	Ace			4.0		100466	Type 44	
	Ace			5.0		100466	Type 44	
SOUTHERN IMPLANTS	External Hex					100353	Type 01	
	Extra-oral			3.75		100353	Type 01	
	IT Connection					100350	Type 04	
	MAX			7.0		100353	Type 01	
	MAX			8.0		100353	Type 01	
	MAX			9.0		100353	Type 01	
	MAX-IT			7.0		100380	Type 07	
	MAX-IT			8.0		100380	Type 07	
	MAX-IT			9.0		100380	Type 07	
	Tri-MAX			7.0		100412	Type 24	
	Tri-MAX			8.0		100462	Type 40	
	Tri-MAX			9.0		100462	Type 40	
	Tri-Nex			3.5		100359	Type 12	
	Tri-Nex			4.3		100361	Type 13	
	Tri-Nex			5.0		100361	Type 13	
STERNGOLD DENTAL	External Hexagon		RP	3.75		100353	Type 01	
	External Hexagon		RP	4.0		100353	Type 01	
	External Hexagon		RP	5.0		100353	Type 01	
	Stern IC Internal Octagon		RN	4.1		100350	Type 04	
STRAUMANN	Bone Level Implants Roxolid SLActive		NC	3.3		100484	Type 53	
	Bone Level Implants Roxolid SLActive		RC	4.1		100485	Type 54	
	Bone Level Implants Roxolid SLActive		RC	4.8		100485	Type 54	
	Bone Level Implants Titanium SLA		NC	3.3		100484	Type 53	
	Bone Level Implants Titanium SLA		RC	4.1		100485	Type 54	
	Bone Level Implants Titanium SLA		RC	4.8		100485	Type 54	
	Bone Level Tapered Implants Roxolid SLA	021.3508	NC	3.3	8	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3510	NC	3.3	10	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3512	NC	3.3	12	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3514	NC	3.3	14	100484	Type 53	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New	
STRAUMANN	Bone Level Tapered Implants Roxolid SLA	021.3516	NC	3.3	16		100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3518	NC	3.3	18		100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.5508	RC	4.1	8		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5510	RC	4.1	10		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5512	RC	4.1	12		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5514	RC	4.1	14		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5516	RC	4.1	16		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5518	RC	4.1	18		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7508	RC	4.8	8		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7510	RC	4.8	10		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7512	RC	4.8	12		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7514	RC	4.8	14		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7516	RC	4.8	16		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7518	RC	4.8	18		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.3308	NC	3.3	8		100484	Type 53	
	Bone Level Tapered Implants Roxolid SLActive	021.3310	NC	3.3	10		100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3312	NC	3.3	12		100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3314	NC	3.3	14		100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3316	NC	3.3	16		100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3318	NC	3.3	18		100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.5308	RC	4.1	8		100485	Type 54	
	Bone Level Tapered Implants Roxolid SLActive	021.5310	RC	4.1	10		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5312	RC	4.1	12		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5314	RC	4.1	14		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5316	RC	4.1	16		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5318	RC	4.1	18		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7308	RC	4.8	8		100485	Type 54	
	Bone Level Tapered Implants Roxolid SLActive	021.7310	RC	4.8	10		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7312	RC	4.8	12		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7314	RC	4.8	14		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7316	RC	4.8	16		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7318	RC	4.8	18		100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive		SC	2.9			100404	Type 22	
	Bone Level Tapered Implants Titanium SLA	021.3408	NC	3.3	8		100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3410	NC	3.3	10		100484	Type 53	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New	
STRAUMANN	Bone Level Tapered Implants Titanium SLA	021.3412	NC	3.3	12		100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3414	NC	3.3	14		100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3416	NC	3.3	16		100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3418	NC	3.3	18		100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.5408	RC	4.1	8		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5410	RC	4.1	10		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5412	RC	4.1	12		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5414	RC	4.1	14		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5416	RC	4.1	16		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5418	RC	4.1	18		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7408	RC	4.8	8		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7410	RC	4.8	10		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7412	RC	4.8	12		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7414	RC	4.8	14		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7416	RC	4.8	16		100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7418	RC	4.8	18		100485	Type 54	New
	Extra Oral						100353	Type 01	
	Narrow Neck						100390	Type 17	
	Ortho						100368	Type 09	
	Palatal						100368	Type 09	
	Standard Implants Roxolid SLActive		RN	3.3			100350	Type 04	
	Standard Implants Roxolid SLActive		RN	4.1			100350	Type 04	
	Standard Implants Roxolid SLActive		RN	4.8			100350	Type 04	
	Standard Implants Roxolid SLActive		WN	4.8			100350	Type 04	
	Standard Implants Titanium SLA		RN	3.3			100350	Type 04	
	Standard Implants Titanium SLA		RN	4.1			100350	Type 04	
	Standard Implants Titanium SLA		RN	4.8			100350	Type 04	
	Standard Implants Titanium SLA		WN	4.8			100350	Type 04	
	Standard Plus Implants Roxolid SLActive		NNC	3.3			100560	Type 68	
	Standard Plus Implants Roxolid SLActive		RN	3.3			100350	Type 04	
	Standard Plus Implants Roxolid SLActive		RN	4.1			100350	Type 04	
	Standard Plus Implants Roxolid SLActive		RN	4.8			100350	Type 04	
	Standard Plus Implants Roxolid SLActive		WN	4.8			100350	Type 04	
	Standard Plus Implants Titanium SLA		RN	3.3			100350	Type 04	
	Standard Plus Implants Titanium SLA		RN	4.1			100350	Type 04	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
STRAUMANN	Standard Plus Implants Titanium SLA		RN		4.8		100350	Type 04
	Standard Plus Implants Titanium SLA		WN		4.8		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		RN		3.3		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		RN		4.1		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		WN		4.8		100350	Type 04
	Tapered Effect Implants Titanium SLA		RN		3.3		100350	Type 04
	Tapered Effect Implants Titanium SLA		RN		4.1		100350	Type 04
	Tapered Effect Implants Titanium SLA		WN		4.8		100350	Type 04
SWEDEN & MARTINA	EZ						100380	Type 07
	Global				3.8		100479	Type 48
	Global				4.3		100479	Type 48
	Global				4.8		100479	Type 48
	Global				5.5		100479	Type 48
	Kohno				3.8		100440	Type 32
	Kohno				4.25		100380	Type 07
	Kohno				5.0		100380	Type 07
	Premium - new				3.3		100440	Type 32
	Premium - new				3.8		100440	Type 32
	Premium				4.25		100380	Type 07
	Premium				5.0		100380	Type 07
	Premium - old				3.3		100402	Type 21
SYBRON	Endopore External				3.5		100436	Type 30
	Endopore External				4.1		100353	Type 01
	Endopore External				5.0		100353	Type 01
	Entegra External				3.25		100353	Type 01
	Entegra External				4.0		100353	Type 01
	Pitt-Easy				4.0		100412	Type 24
	Pitt-Easy				4.9		100412	Type 24
	SybronPRO XRT				3.3		100378	Type 06
	SybronPRO XRT				3.5		100378	Type 06
SybronPRO XRT				4.1		100378	Type 06	
TA-DENT					3.3		100455	Type 38

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
TA-DENT				4.8		100378	Type 06	
				5.6		100480	Type 49	
TAG DENTAL	Axis			3.3		100431	Type 27	New
	Axis			3.75		100431	Type 27	New
	Axis			4.2		100431	Type 27	New
	Axis			5.0		100431	Type 27	New
	Axis			6.0		100431	Type 27	New
	Massif			3.75		100431	Type 27	New
	Massif			4.2		100431	Type 27	New
	Massif			5.0		100431	Type 27	New
	Massif			6.0		100431	Type 27	New
	TBR	External Hexagon					100353	Type 01
Internal Octagon						100398	Type 18	
M-implant				3.2		100402	Type 21	
M-implant				3.9		100402	Type 21	
M-implant				4.7		100402	Type 21	
Z1-M implant				3.2		100402	Type 21	
Z1-M implant				3.9		100402	Type 21	
Z1-M implant				4.7		100402	Type 21	
TEKKA	Inkone			3.5		100480	Type 49	
	Inkone			4.5		100480	Type 49	
	Krestal			3.5		100378	Type 06	
	Krestal			4.0		100378	Type 06	
	Krestal			4.2		100378	Type 06	
	Krestal			5.1		100378	Type 06	
	Progress			3.4		100378	Type 06	
	Progress			4.0		100378	Type 06	
	T-quest			4.0		100480	Type 49	
	T-quest			4.5		100480	Type 49	
	Twinkone			3.5		100390	Type 17	
TFI SYSTEM	TFI			3.3		100425	Type 26	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New	
TFI SYSTEM	TFI			3.75			100425	Type 26	
	TFI			4.2			100425	Type 26	
	TFI			5.0			100425	Type 26	
THOMMEN MEDICAL	SPI System Platform			3.5			100404	Type 22	
	SPI System Platform			4.0			100402	Type 21	
	SPI System Platform			4.5			100386	Type 15	
	SPI System Platform			5.0			100386	Type 15	
	SPI System Platform			6.0			100386	Type 15	
TIDAL SPIRAL				3.75			100431	Type 27	
				4.2			100431	Type 27	
				5.0			100431	Type 27	
				6.0			100431	Type 27	
TRATE	Rootform			4.0			100497	Type 57	New
TRI DENTAL IMPLANTS	TRI-Octa 4.8			3.75			100544	Type 61	
	TRI-Octa 4.8			4.1			100544	Type 61	
	TRI-Octa 4.8			4.7			100544	Type 61	
	TRI-Vent			3.5			100440	Type 32	
U-IMPL	External Hex						100353	Type 01	
	Internal Hex						100440	Type 32	
WARANTEC	Oneplant Implant			4.3			100412	Type 24	
	Oneplant Internal Tapered						100378	Type 06	
WEGO				3.4	9		100484	Type 53	
				4.3	7		100455	Type 38	
				4.3	9		100411	Type 23	
WIELAND IMPLANT SYSTEM				3.5			100472	Type 45	
				4.3			100472	Type 45	
				5.0			100472	Type 45	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	SmartPeg	Art. No	New
YAMAHACHI DENTAL	μ-one HA Implant 2 piece			3.4		100455	Type 38	
	μ-one HA Implant 2 piece			3.8		100455	Type 38	
	μ-one HA Implant 2 piece			4.3		100431	Type 27	
	μ-one HA Implant 2 piece			5.0		100398	Type 18	
ZIMMER DENTAL	AdVent			4.5		100446	Type 35	
	AdVent			5.7		100444	Type 34	
	Eztetic Implant System CT31			3.1		100543	Type 60	
	Screw-vent		3.5	3.3		100440	Type 32	
	Screw-vent		3.5	3.7		100440	Type 32	
	Screw-vent		4.5	4.7		100431	Type 27	
	Spline			3.25		100436	Type 30	
	Spline			3.75		100436	Type 30	
	Spline			5.0		100436	Type 30	
	SwissPlus			3.8		100425	Type 26	
	SwissPlus			4.8		100350	Type 04	
	Tapered Screw-vent		3.5	3.7		100440	Type 32	
	Tapered Screw-vent		3.5	4.1		100440	Type 32	
	Tapered Screw-vent		4.5	4.7		100431	Type 27	
Tapered Screw-vent		5.7	6.0		100442	Type 33		
ZUGA MEDICAL	Bone Level			3.5		100554	Type 63	
	Bone Level			5.0		100557	Type 65	

Implant Company	Abutment Model	Catalog No	Platform	Angle	Diameter (mm)	Length or GH (mm)	SmartPeg	Art. No	New
CAMLOG	Comfour Bar abutment		Yellow		3.8		100579	Type 73	New
	Comfour Bar abutment		Purple		4.3		100579	Type 73	New
	Comfour Bar abutment		Blue		5.0		100580	Type 74	New
COCHLEAR	Baha BA 210					5.5	100487	Type 55	
	Baha BA 210					8.5	100487	Type 55	
	Baha BA 300, BIA 300					6.0	100487	Type 55	
	Baha BA 300, BIA 300					9.0	100487	Type 55	
	Baha BA 400, BIA 400					10.0	100487	Type 55	
	Baha BA 400, BIA 400					12.0	100487	Type 55	
	Baha BA 400, BIA 400					14.0	100487	Type 55	
	Baha BA 400, BIA 400					6.0	100487	Type 55	
	Baha BA 400, BIA 400					8.0	100487	Type 55	
	Vistafix Standard abutment					3.0	100363	Type A3	
	Vistafix Standard abutment					4.0	100363	Type A3	
	Vistafix Standard abutment					5.5	100363	Type A3	
	Vistafix Standard abutment					7.0	100363	Type A3	
	Vistafix VXA 300					3.5	100363	Type A3	
	Vistafix VXA 300					4.5	100363	Type A3	
	Vistafix VXA 300					6.0	100363	Type A3	
	Vistafix VXA 300					7.5	100363	Type A3	
NEOSS	Access abutment						100483	Type 52	
NOBEL BIOCARE	Brånemark system MUA		NP				100363	Type A3	
	Brånemark system MUA		RP				100363	Type A3	
	On1 Base CC		NP				100579	Type 73	New
	On1 Base CC		RP				100579	Type 73	New
	On1 Base CC		WP				100579	Type 73	New
	Mua		NP				100415	Type 25	New
	Mua		RP				100415	Type 25	New
	Mua		WP				100415	Type 25	New
OTICON MEDICAL	Ponto Implant abutment					4.0	100487	Type 55	

Implant Company	Abutment Model	Catalog No	Platform	Angle	Diameter (mm)	Length or GH (mm)	SmartPeg	Art. No	New	
STRAUMANN	Coping for screw-retained abutment, Ti Bridge	023.2749	NC		3.5			100415	Type 25	New
	Coping for screw-retained abutment, Ti Bridge	023.4751	NC/RC		4.6			100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2749	NC	17	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2755	NC	17	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0010	NC	17	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2751	NC	30	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2757	NC	30	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0012	NC	30	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4747	RC	17	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4752	RC	17	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4749	RC	30	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4754	RC	30	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0014	RC	17	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0016	RC	30	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2750	NC	17	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2756	NC	17	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0011	NC	17	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2752	NC	30	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2758	NC	30	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0013	NC	30	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4748	RC	17	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4753	RC	17	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4750	RC	30	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4755	RC	30	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0017	RC	30	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0015	RC	17	4.6	5.5		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2745	NC	0	3.5	1		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2746	NC	0	3.5	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2753	NC	0	3.5	4		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2747	NC	0	4.6	1		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2748	NC	0	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2754	NC	0	4.6	4		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.4745	RC	0	4.6	1		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.4746	RC	0	4.6	2.5		100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.4751	RC	0	4.6	4		100415	Type 25	New

Implant Company	Abutment Model	Catalog No	Platform	Angle	Diameter (mm)	Length or GH (mm)	SmartPeg	Art. No	New
THOMMEN MEDICAL	SPI System Bar abutment						100447	Type 36	
	VARIOmulti abutment						100415	Type 25	



A fresh start every time!

To avoid damage to the implant, the SmartPeg material is quite soft, which limits the lifetime of the threads. Therefore, SmartPegs are single use items.

www.osstell.com/webshop

