




Zuordnung Spannzangen der W&H Labor-Handstücke

Correlation chucks of the W&H laboratory handpieces

MF-PERFECTA MF-TECTORQUE




REF 03201500
Spannzange Ø 2,35 mm
Chuck Ø 2.35 mm

REF 03202600
Spannzange Ø 3,00 mm
Chuck Ø 3.00 mm

MF-PERFECTA MF-TECTORQUE



perfecta 300/600
perfecta 900



REF 03205500
Spannzange Ø 2,35 mm
Chuck Ø 2.35 mm

REF 03205600
Spannzange Ø 3,00 mm
Chuck Ø 3.00 mm