



The formula for a
long instrument life



f1 service
oil

The More Intelligent Formula

Prolong the lifespan of your instruments

For Every Application

W&H Service Oil F1 is a product of W&H research. It has been developed hand in hand with modern W&H instrument technology. This oil gives you the confidence of knowing that you are providing your W&H instruments with the best possible care:

- > Turbines
- > Straight and contra-angle handpieces
- > Surgical instruments
- > Air scalers
- > Air motors



Perfect for manual maintenance



W&H Service Oil F1 in the MD-400 spray can



W&H Service Oil F1 in the MD-500 bottle for use in the Assistina

W&H dental drives achieve top speeds of up to 390,000 r.p.m. F1 provides excellent lubrication

and internal cleaning thereby ensuring fault-free operation of the instruments.

With the test champion oil for all standard drives

The More Intelligent Service Oil

W&H Service Oil F1 is an all synthetic, special oil of maximum purity, which has been specially developed for the demanding and extreme requirements of dental instruments. Compared with other well-known brand oils, the Service Oil F1 was subjected to a series of rigorous tests.

The Result: Test winner in all significant performance features. The test result: The better the oil, the lower the wear. (Indentation diameter of friction WK [mm] determined from TU-Wien according to the SRV-test [translatory oscillation-test] standard DIN 51 834.)

Sterilization resistant



F1 is of course fully sterilizable, absolutely resistant to hydrolysis, stable against oxidation and is shown to be highly compatible with disinfectants.

Physiologically safe

F1 contains no toxic constituents, is physiologically harmless and environmentally compatible. The substances it contains comply with the applicable international regulations.


High material compatibility

F1 also safeguards the durability of plastic parts, coatings, seals, glass etc.

Enormously economical

Your instruments will thank you for regular maintenance with F1 by requiring fewer repairs. The advantages: considerable savings in cost and increased profitability of your W&H instruments.

Short lubrication intervals – long service life



Intensive maintenance will considerably increase the life of your instruments.

Intelligent oil + short maintenance intervals = longer lifespan

Diagrams, symbol photos, additional equipment and the content of illustrated accessories do not form part of the scope of supply.



Oil spray can

Type MD-400

Contents 400 ml

also available in packs of 6



Bottle for Assistina

Type MD-500

Contents 500 ml

Order your Service Oil F1 through your dental supplier.
Suitable spray caps can be found on the Internet at wh.com

Manufacturer:

W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, Austria
t + 43 6274 6236-0
f + 43 6274 6236-55
office@wh.com
wh.com

Distribution:

W&H (UK) Limited
6 Stroud Wood Business Centre
Park Street, St Albans, Hertfordshire AL2 2NJ, United Kingdom
t +44 1727 874990
f +44 1727 874628
office.uk@wh.com
wh.com

