



# LED upgrade set



#### Advantages:

- > Daylight quality LED light
- > Simple installation even with existing units
- > Ten-times longer lifespan than halogen bulbs

Change your existing halogen motor over to daylight quality LED light. The W&H LED upgrade set allows you to upgrade the W&H EA-40 LT motor with just two new components. So you can enjoy daylight quality light generated by innovative LED technology, while still using conventional glass rod instruments.

led

Diagrams, symbol photos, additional equipment and the content of illustrated accessories are not included in delivery.

## Connection diagram



## Technical Data

### LED

for installation in the motor

Colour temperature: 5,500 K

Sterilizable: up to 135°C

### LED electronics system

for installation in the dental unit

Dimensions:

28 mm x 20 mm x 25 mm

Upgrade your EA-40 LT motor from a halogen light to an LED. You will be amazed at just how much more you can see.

No matter what glass rod instruments you use, a motor which has been upgraded to an LED gives you perfect white daylight quality light throughout the entire treatment area.

LED technology has a significantly longer lifespan than conventional halogen bulbs.

What's more, it provides an almost unimaginable intensity of daylight quality light, which is many times brighter than standard halogen light.



The circuit board has been safety tested for its compliance with standard IEC 60601-1 and carries the UL quality mark. It is the responsibility of the specialist customer to ensure that the entire system complies with all applicable standards and regulations.

**W&H Dentalwerk Bürmoos GmbH**

Ignaz-Glaser-Straße 53, Postfach 1

5111 Bürmoos, Austria

t + 43 6274 6236-0

f + 43 6274 6236-55

office@wh.com

wh.com



ISO 13485  
93/42/EEC - Annex II