

W&H Couplings

Versatile system interfaces



Flexible – for use with all standardised turbine tubing. **Full compatibility** – compatible with all turbines, air motors and scalers with Roto Quick and Multiflex®* connections. **360° instrument rotation** – for easy handling, prevents the tube from twisting. **Continuous regulation of the spray water supply** – quick, customised adjustment by the user even –during treatment. **Non-retraction valve** – avoids back-flow of contaminated spray water to minimise the risk of infection.

Diagrams, symbol photos, additional equipment and the content of illustrated accessories are not included in delivery.



Roto Quick couplings:	RQ-24	RQ-34	RQ-14	RQ-03	RQ-04
Light:	yes	yes	–	–	–
Spray control:	–	yes	yes	–	–
Connection:	6-hole	6-hole	Standard 4-hole	2 (3)-hole	Standard 4-hole

W&H Roto Quick couplings

Light, flexible and sterilizable: W&H Roto Quick couplings are versatile to use and are suitable both for instruments with LED+ and LED, as well as halogen light. 24 months of warranty.

W&H coupling for Multiflex®* connection:	RM-34 LED
Light:	LED
Spray control:	yes
Connection:	6-hole



W&H coupling RM-34 LED for Multiflex®* connection

Simple and uncompromising: The W&H coupling RM-34 LED is ready for immediate use. This means: The coupling with LED can be used alongside the conventional halogen couplings without requiring any changes at all in the dental unit. Long service life – 24 months of warranty.

* Multiflex® is a registered trademark of Kaltenbach & Voigt GmbH, Germany

Manufacturer:

W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, Austria
t +43 6274 6236-0
f +43 6274 6236-55
office@wh.com
wh.com

Distribution:

W&H (UK) Limited
6 Stroud Wood Business Centre
Park Street, St Albans, Hertfordshire AL2 2NJ, United Kingdom
t +44 1727 874990
f +44 1727 874628
office.uk@wh.com
wh.com

