



Compact. Lightweight. led.



NEW

Air motors

AM-25

Technical data



Scan the QR code with your mobile phone for more information on air motors

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.



Type	AM-25 L RM	AM-25 A RM AM-25 E RM AM-25 RM	AM-25 A BC AM-25 E BC AM-25 BC
Coupling hose end:	Fixed connection Standard 6-hole	Fixed connection Standard 4-hole	Fixed connection Borden 2(3)-hole
Light (LED version):	LED	–	–
Spray:	internal	internal external –	internal external –
Weight:	77 g	78 g 75 g 73 g	78 g 74 g 73 g
Power:		20 W	
Speed range (min ⁻¹) at an operating pressure of:		2.5 bar: 5,000 – 20,000* 3 bar: 5,000 – 25,000*	
Direction of rotation:		forward / reverse operation mode	

* The power and speed of air drives depend on the quality of the turbine hose and may differ from the specified value.




LED light

The LED in the air motor guarantees optimal view and illumination of the treatment site.



Excellent ergonomics

Easy 360° rotation of the contra-angle handpiece on the motor and the short, lightweight design ensure fatigue-free operation.



Fast removal of tooth structure

Extremely powerful motor for efficient operation with high torque.

W&H Dentalwerk Bürmoos GmbH
 Ignaz-Glaser-Straße 53, Postfach 1
 5111 Bürmoos, Austria
 t + 43 6274 6236-0
 f + 43 6274 6236-55
 office@wh.com
 wh.com