



# Power Demand Piezo-Tips

| Tip         |           | Piezo-Scaler |         |         |         |                |         |         |         |         |         |
|-------------|-----------|--------------|---------|---------|---------|----------------|---------|---------|---------|---------|---------|
|             |           | Pyon1        |         | Pyon2   |         | Piezo-Built in |         | Tigon   |         | Tigon+  |         |
|             |           | Optimum      | Maximum | Optimum | Maximum | Optimum        | Maximum | Optimum | Maximum | Optimum | Maximum |
| Scaling     | 1U        | 30           | 40      | 30      | 40      | 30             | 40      | 30      | 40      | 30      | 40      |
|             | 2U        | 20           | 30      | 20      | 30      | 20             | 30      | 20      | 30      | 20      | 30      |
|             | 3U        | 20           | 30      | 20      | 30      | 20             | 30      | 20      | 30      | 20      | 30      |
|             | 4U        | *            | *       | 30      | 40      | 30             | 40      | 30      | 40      | 30      | 40      |
| Perio       | 1P        | 10           | 30      | 10      | 30      | 10             | 30      | 10      | 30      | 10      | 30      |
|             | 2Pr / 2PI | *            | *       | 20      | 30      | 20             | 30      | 20      | 30      | 20      | 30      |
|             | 3Pr / 3PI | *            | *       | *       | *       | *              | *       | 20      | 30      | 20      | 30      |
| Implant     | 1I        | *            | *       | 25      | 30      | 25             | 30      | 25      | 30      | 25      | 30      |
| Endo        | 1R        | *            | *       | 20      | 40      | 20             | 40      | 20      | 40      | 20      | 40      |
|             | 1E        | *            | *       | 15      | 15      | 15             | 15      | 15      | 15      | 15      | 15      |
|             | 2E        | *            | *       | 35      | 35      | 35             | 35      | 35      | 35      | 35      | 35      |
|             | 3E        | *            | *       | 20      | 20      | 20             | 20      | 20      | 20      | 20      | 20      |
| Restoration | 1R        | *            | *       | 35      | 40      | 35             | 40      | 35      | 40      | 35      | 40      |
|             | 2R        | *            | *       | 35      | 40      | 35             | 40      | 35      | 40      | 35      | 40      |
|             | 3Rd / 3Rm | *            | *       | 35      | 40      | 35             | 40      | 35      | 40      | 35      | 40      |

\* not applicable